




GEMINI 2D

INTERFEROMETER



(주)원우시스템즈 - NIROES 한국공식대리점
07072 서울시 동작구 신대방 1가길 38, 106동 209호
TEL : 02) 533-6720, FAX : 02)3289-1293
<https://www.wonwoosystem.co.kr>

GEMINI 2D *for time-resolved measurements*



The GEMINI 2D is the advanced model of the GEMINI interferometer. It is a compact and ultra-stable interferometer, capable of generating a pair of **collinear and phase-locked replicas of ultra-short pulses**, with unrivaled stability and robustness.

The GEMINI 2D is specially designed to:

- keep constant the **dispersion** introduced during the scan of the delay between the two replicas
 - keep fixed the **absolute arrival time** of one of the two replicas (with attosecond stability)
- High throughput
- 1 attosecond stability
- Scan range selectable by the user
- Compact and insensitive to vibrations


Main Applications

The GEMINI 2D is the ideal device for **time-resolved measurements**, where it is crucial to preserve the pulse duration and the synchronization with other light pulses.


Perfectly suitable for **two-dimensional (2D) spectroscopy experiments**, to generate a pair of collinear and phase-locked pump pulses



Two-dimensional Spectroscopy (in pump-probe geometry)




Bidimensional maps of 2D Electronic Spectroscopy measurements obtained on a Light Harvesting (LH1) complex of a sample of Rhodospirillum Rubrum for three different population times T^{\wedge} (15 fs, 45 fs, 4000 fs).



- The GEMINI 2D interferometer is placed in the pump beam before the sample, allowing one to generate two collinear and phase-locked pump pulses (1 and 2) with a relative variable delay T^{\wedge} (coherence time).
- The absolute arrival time of pulse 2 is kept fixed in order not to change the delay T^{\wedge} (population time) during the scan of T^{\wedge} .
- In the pump-probe geometry, the rephasing and non-rephasing signals - P_{NR+R} - are emitted in the same direction as the probe (3), which also acts as a local oscillator (LO) that heterodynes them on the detector.


Performances

REPEATABILITY TEST



Interference patterns between the two replicas generated by the GEMINI 2D Interferometer, measured for different scan number indicated on the right. The lower one is the average of 40 scans. The measurements show an **exceptional repeatability!**

STABILITY TEST



(a) Fringe pattern generated by the interference of the two replicas created by the GEMINI 2D Interferometer; (b) sequence of fringe patterns acquired by keeping the relative delay fixed for 30 min; and (c) corresponding phase fluctuations. The **phase jitter is less than 1/360th of the optical cycle.**

Technical Specifications

VERSION		S	L
Spectral range [nm]		400 - 2300 (Standard) 250 - 3500 (Ultra-broadband)	
Max. Delay τ [fs @ $\lambda=600$ nm]	SYM	± 400	± 1050
	ASY	-100 \rightarrow 700	-100 \rightarrow 2000
Delay τ Stability		< 1 attosecond	
Dimensions [mm]		180 x 180 x 90	
Weight [kg]		2	

Spectral Resolution

