Integrating Spheres

Integrating spheres are unique tools used in many laser & photonics applications. Gigahertz-Optik offers a wide range of integrating spheres for both detection purposes and for uniform light sources. To meet the many varied application demands, spheres with different diameters are offered as well as different coatings and measurement geometries. Besides the modular UM series for light measurements and uniform light source set-ups Gigahertz-Optik offers the UP series with precisely aligned port geometries for reflectance and transmittance measurements. Complete systems can be assembled using accessory components and an instrument from Gigahertz-Optik's photometer, radiometer and colorimeter program. To satisfy specific application demands, custom design and private labeling solutions are offered.

Table of Contents						Page
UM Series	Descriptio	n and exam	ples for inc	lividual UM series	integrating sphere set-ups	95-99
Model		Sphere		Accessory	Description	Specifi-
	Туре	Diameter mm / in	Coating	Туре		cations
UMBB-100	Spun Al	100 / 4	BaSO4		Basic integrating spheres for detection & light source set-ups.	100
UMBB-150	Spun Al	150 / 6	BaSO4		First, port frames, baffles and stands are added to the spun	
UMBB-210	Spun Al	210/8.3	BaSO4		Then light sources variable apertures detectors sample hold-	
UMBB-300	Spun Al	300 / 12	BaSO4		ers, auxiliary lamps and other components can be assembled	
UMBB-500	Spun Al	500 / 20	BaSO4		to the port frames per end-user specifications. Finally, periph-	
UMBB-1000	Spun Al	1000 / 40	BaSO4		eral instruments like optometers and lamp power supplies com-	
UMBB-1700	Spun Al	1700/67	BaSO4		Three different coatings are supplied: barium sulfate (UMBB	
UMBG-100	Spun Al	100 / 4	Gold		models) can be used for all sizes and most accessories, gold	
UMBG-150	Spun Al	150 / 6	Gold		OMBG models) for Infrared applications and OP.DI.MA. (UMBK models) for UV-VISIBLE-IB broadband applications with limited	
UMBG-300	Spun Al	300 / 12	Gold		availability sphere sizes and accessories.	
UMBK-190	Spun Al	190 / 7.5	OPDIMA		Please contact the factory or your local Gigahertz-Optik repre-	
UMBK-250	Spun Al	250 / 10	OPDIMA		sentative for configuration assistance.	
UMBK-460	Spun Al	460 / 18	OPDIMA			
UMSS-SM14/M6				Sphere Mount	Post mount socket for small diameter spheres	101
UMSS-SMT				Sphere Mount	Flange mounted adapter for UMSS-BT use	101
UMSS-SMV				Sphere Mount	Flange mounted adapter for post mount use	101
UMSS-BT150				Sphere Stand	Bench top stand for UMBB-150 & UMBG-150	101
UMSS-BT210				Sphere Stand	Bench top stand for UMBK-190 & UMBB-210	101
UMSS-BT300				Sphere Stand	Bench top stand for UMBB-300, UMBG-300, UMBK-250	101
UMSS-BT500				Sphere Stand	Bench top stand for UMBK-460 & UMBB-500	101
UMSS-BT1000				Sphere Stand	Bench top stand for UMBB-1000	101
UMSS-BT1700				Sphere Stand	Bench top stand for UMBB-1700	101
UMSS-HF300				Sphere Stand	Hinge frame bench top stand for UMBB-300	102
UMSS-HF500				Sphere Stand	Hinge frame bench top stand for UMBB-500	102
UMSS-HF1000				Sphere Stand	Hinge frame stand for UMBB-1000	102
UMSS-HF1700				Sphere Stand	Hinge frame stand for UMBB-1700	102
UMPF-xx				Port Frames	From 0.5 to 5.0 inch dia Fix mounted to UMB	101
UMPP-xx				Port Plugs	From 0.5 to 5.0 inch to plug in UMPF-0.5 to 5.0	103
UMPR-xx				Port Reducer	From 0.5 to 5.0 inch	103
UMPD-37				Detector Port	For 37 mm dia. light detectors. Fix mounted to UMBB	104
UMPD-45				Detector Port	For 45 mm dia. light detectors. Fix mounted to UMBB	104
UMPF-LSOK30				Port Frame	For LS-OK30 light source assembling	104
UMPF-LSOK45				Port Frame	For LS-OK45 light source assembling	104
UMPB-SM				Port Baffle	For light baffling between ports. Sphere wall mounted	104
UMPB-PM				Port Baffle	For light baffling to sphere center. Post mounted	104
UMPB-IL				Port Baffle	For light baffling of in-line ports. Spider mounted	104
UMPA-0.5/11				Port Adapter	To assemble 11-type detectors to UMPF-0.5	105-106
UMPA-0.5/11RD				Diffuser	RADIN diffuser to be assembled in UMPA-0.5/11	105-106
UMPA-0.5/F				Port Adapter	Universal machineable. Fits to UMPF-0.5	105-106
UMPA-1.0/F				Port Adapter	Universal machineable. Fits to UMPF-1.0	105-106

Page 94 - Integrating Spheres



Table of Contents

UM Series					Page	
Model		Sphere		Accessory	Description	Specifi-
	Туре	Diameter mm / in	Coating	Туре		cations
UMPA-1.5/LGA13				Port Adapter	Light guide adapter for Spectrometer. Fits UMPF-1.5 (1.5")	105-106
UMPA-2.0/AXL				Port Adapter	For auxiliary lamp set-up with LS-IS1.5. Fits UMPF-2.0 (2.0")	105-106
UMPA-2.0/LED5				Port Adapter	For LED-5xx LED test sockets. Fits UMPF-2.0 (2.0")	105-106
UMPA-4.0/QD1.9				Port Adapter	48 mm dia. diffuse quartz window. Fits UMPF-4.0 (4.0")	105-106
UMPA-4.0/QD3.0				Port Adapter	75 mm dia. diffuse quartz window. Fits UMPF-4.0 (4.0")	105-106
UMPA-5.0/LED5				Port Adapter	For LED-5xx LED test sockets. Fits UMPF-5.0 (5.0")	105-106
UMPA-5.0/SH				Port Adapter	Sample base holder with 4" open aperture. Fits UMPF-5.0 (5.0")	105-106
UMPA-5.0/SH-P				Port Adapter	Port plug for UMPA-5.0/SH	105-106
UMPA-5.0/SH-CM				Port Adapter	Port plug for UMPA-5.0/SH with universal sample mount	105-106
UMLA-1.5B				Lamp Adapter	Base for center position lamp adapter	106
UMLA-300				Lamp Adapter	For center lamp/sample position inside UMBB-300	106
UMLA-500				Lamp Adapter	For center lamp/sample position inside UMBB-500	106
UP Series						107
UPB-50-L	CNC AI	50/2	BaSO4		Precision machined single port sphere	108
UPB-100-F	CNC AI	100 / 4	BaSO4		Precision machined multi-port flexible design sphere with plugs	108
UPB-150-F	CNC AI	150 / 6	BaSO4		Precision machined multi-port flexible design sphere with plugs	109
UPB-250-F	Spin Al	250 / 10	BaSO4		Precision machined multi-port flexible design sphere with plugs	109
UPB-500-F	Spin Al	500 / 20	BaSO4		Precision machined multi-port flexible design sphere with plugs	110
UPK-30-L	CNC AI	30 / 1.2	OPDIMA		Precision machined single port sphere	111
UPK-50-F	CNC AI	50/2	OPDIMA		Precision machined multi-port flexible design sphere with plugs	111-112
UPK-50-L	CNC AI	50/2	OPDIMA		Precision machined single port sphere	111-112
UPK-50-ULS	CNC AI	50/2	OPDIMA		Precision machined (8) LED port sphere. Uniform light source.	111-112
UPK-100-F	CNC AI	100 / 4	OPDIMA		Precision machined multi-port flexible design sphere with plugs	112
UPK-100-L	CNC AI	100 / 4	OPDIMA		Precision machined single port sphere	112
UPK-150-F	CNC AI	150 / 6	OPDIMA		Precision machined multi-port flexible design sphere with plugs	113
UPK-190	CNC AI	190 / 7.5	OPDIMA		Precision machined sphere. Custom design configuration.	113
UPK-250	CNC AI	250/10	OPDIMA		Precision machined sphere. Custom design configuration.	114
UPG-75-G	CNC AI	75/3	Gold		Precision machined multi-port flexible design sphere with plugs	114
UPG-100-F	CNC AI	100 / 4	Gold		Precision machined multi-port flexible design sphere with plugs	115
Accessories						
ULT-18				Light Trap	For assembly to 11-type ports	115
UFC				Fiber Socket	FC, SC, SMA and ST fiber coupler adapter for 11-type ports	116
Light Sources & Li	ght Sourc	e Accesso	ries			
LS-IS1.5				Internal Source	Lamp housing for 5 to 100 W QH-lamps. Fits UMPF-1.5 (1.5")	117
LS-OK30				External Source	Lamp housing for 5 to 100 W QH-lamps. Fits UMPF-LSOK30	118
LS-OK30-FH				Filter Holder	For use with LS-OK30	118
LS-OK30-FHI				Filter Holder	Insert for 2 inch diameter filter for use with LS-OK30-FH	118
LS-OK30-MIR				Variable Aperture	High resolution Iris diagraph for use with LS-OK30	118
LS-OK30-VA				Variable Aperture	for use with LS-OK30	118
LS-OK45				External Source	Lamp housing for 250 W QH-lamps.	119
LS-OK45-LG				Light Guide	Adapter. Color temperature control operation with CT-3701	119
VAM-45				Variable Aperture	Stepper motor operated variable aperture. Needs SMC-01.	119
FWM-12/2				Filter Wheel	Stepper motor controlled filter wheel. Needs SMC-01.	119
LS-CB				Beam Source	Collimated light source for reflectance & transmittance meas.	119
LPS-250				Power Supply	Precision lamp power supply (up to 250 W lamps)	120
LCRT-2000				Power Supply	Lamp power supply (up to 10 W lamps)	120
SMC-01				Motor Controller	4-channel stepper motor controller with interface.	121
Application Service	9					
Gigahertz-Optik offe	rs full appli	cation servi	ce for integ	rating spheres and	d integrating sphere systems.	122

Gigahertz-Optik

UM Series

Years of experience in manufacturing integrating spheres has taught us that a line of standard spheres will not fulfill all possible customer needs. Most often either the end-user accepts a standard sphere as a compromise solution or must have the stock item custom modified adding extra cost and lead time. The new **UM** integrating sphere series allows individual integrating sphere set-up in the specific configuration required with only

those components needed. An added benefit to this design freedom is that you don't pay for something you don't need. With standardized stock components and an innovative assembly process the UM series offers economical prices and short delivery times for integrating spheres that others can only claim as custom products. Port frames with free apertures in different diameters as well as detector and light source ports can be added and positioned on the basic spheres according to end-user design. **UMB Basic Spheres** are available in different diameters and with different coatings.

Sphere Accessory Components such as port plugs, reducers, adapters, baffles, sphere mounts and stands are available. Application Accessory Components & Instruments like fiber optic adapters, light sources & detectors, photometers, radi-

Individual Sphere Design

ometers, colorimeters & calibration standards allow complete systems to be assembled.

The following pages describe typical application set-ups as well as individual component specifications. Light measurement instrumentation and calibration equipment technical specifications and descriptions are shown in other chapters of this catalog.

Sphere Components

Light Detector



UM Sphere Configuration



Assembly Orientation

UMB Basic Spheres are assembled using two hollow hemispheres made from aluminum sheet stock.

The **Orbit Flange** located on open side of each hemisphere is used to mechanically mate the sphere halves. The standard position for the orbit flange is vertical as shown in the above drawing. But its actual position can be specified at any angle from vertical. Note that port frames, sphere stands or other accessory may limit the choice of angle. 5 deg. angle steps with a tolerance of +/- 5 deg. are recommended.

To specify position coordinates for **Accessory Components**, the same 5 deg. +/-5 deg. resolution is recommended. Position can be mapped using the sideview 'SV' and top-view "TV' angle charts as shown above.

The sphere Al sheet forming process only allows a typical position tolerance of +/- 5 degrees.

For more precise alignment **UP** series spheres are recommended.

Page 96 - Integrating Spheres



UM Series "Set-up Examples"

Luminous Flux & Radiant Power Measurement of Spot Sources

The integrating sphere is the preferred tool for the measurement of radiant power or luminous flux of spot sources with a diverging light beam. Typical examples of spot sources are LEDs, LED panels, laser diodes, fiber bundle light sources, reflector spot lamps, flat panel displays, endoscopes and others.

The reason why the integrating sphere is preferable to a planar detector is that typical photodiode detectors are limited in active area, size and acceptance angle. Integrating spheres offer large area measurement apertures and large acceptance angles in observance to the rule that the 'measurement aperture area should not exceed 5% of the total sphere area'. Temperature is another reason for using integrating spheres in applications involving high power tungsten or xenon spot lamps. Large diameter spheres are required to maintain temperature within operating limits, prohibit heat damage and for stability/ accuracy of the measurement.

In this set-up a measurement port and a detector port is needed. Photometric, radiometric or colorimetric detectors can be mounted onto the integrating sphere with an accompanying read-out instrument. In applications where wide variations in sphere housing temperature may occur, detectors with temperature stabilization are recommended to avoid drift of the measured signal. A spectrometer can also be used with or as an alternative to the integral light detector. Baffles are placed between the lamp and detector to prohibit direct illumination of the detector by the incoming signal.

The complete sphere system can be calibrated in absolute quantities such as luminous flux

1

Directional LED Measurement

5

3



(Im), radiant power (W) or spectral radiant power (W/nm). Lamp standards are normally used for calibrations in illuminance (Ix) or spectral irradiance (W/(m² nm)). See *Calibration Standards* sec-



Spot Lamp Measurement



150 mm (6 in) dia sphere with photometric detector for luminous flux measurements.

- 1. UMBB-150: 150 mm (6") dia integrating sphere
- 2. UMSS-SM14: sphere socket
- 3. UMPF-1.0: 1.0" dia port frame
- 4. UMPB-SM: port baffle
- 5. UMPF-0.5: 0.5" dia port frame
- UMPA-0.5/11: port adapter for 11-type light detector
- 7. UMPA-0.5/RADIN: diffuser for UMPA-0.5/11 port adapter
- 8. VL-1101-2: photometric detector with calibration data connector for P-9710-1
- 9. P-9710-1: single channel optometer with RS232 Inter-

face 10.KDW: Calibration of luminous flux sensitivity

Options:

2

- TD-11VL01-2: temp. stabilized photometric detector
- P-2000: two channel optometer with RS232 and IEEE488 interface.
- UMPF-2.0: 2.0" dia port frame
- UMPR-2.0/xx: port reducer plug, xx = aperture dia
- PM-B + PM-P: post stand
- UMSS-BT150: bench-top stand alternative to UMSS-

500 mm (20 in) dia sphere with colorimetric detector for luminous flux and color temperature measurements.

- 1. UMBB-500: 500 mm (20") dia integrating sphere
- 2. UMŠS-SMT
- 3. UMSS-BT500: bench-top stand
- 4. UMPF-4.0: 4.0" dia port frame
- 5. UMPR-4.0/xx: 1 or more port reducer, xx = aperture dia
- 6. UMPB-SM: port baffle
- 7. UMDP-37: detector port for 37type detectors
- 8. CT-3701: luminous flux, color temperature and xy chromaticity value detector. Tempera-

ture stabilized

- 9. P-9801: 8-channel optometer with RS232 and IEEE488 interface
- 10.OS-P9801: remote control software
- 11.KDW: calibration of luminous flux sensitivity and color temperature
- Options:
- UMSS-HF500: hinge Frame bench top stand for full opening of the front hemisphere
- Additional radiometric detectors

UM Series "Set-up Examples"

Luminous Flux & Radiant Power Lamp Measurement

the detector by the incoming signal.

The complete integrating sphere light detection system can be calibrated in absolute optical

measurement quantities such us luminous flux (Im), radiant power (W) or spectral radiant power (W/ nm). Please refer to the *Calibration Standards* section for further



Large-size Lamp Measurement

(88)

5

Measuring total radiant power or luminous flux emitted by lamps is one of the most common applications for integrating spheres. Other than spot sources, lamps emit light in all directions. The integrating sphere collects all of the emitted light entering it and supplies an integrated signal to its mounted light detector.

Lamps come in many different form, size and output power today. Examples include tungsten filament lamps, fluorescent tube & bulb shape lamps, high power xenon sources and many more. Selecting the size of the integrating sphere is typically based on the maximum dimension of the lamp to be measured. The sphere diameter should be at least ten times the maximum dimension of the lamp. For tube lamps it should be twice the largest dimension of the source. For flux measurements a multidirectional lamp is placed in the center of the sphere. Different forms and sizes of sockets and individual lamp holders are required for each different lamp type. For smaller size lamps a top-load port can be used for insertion and removal. For large size lamps integrating spheres with a hinged frame allowing one half of the sphere to fully open are recommended.

Photometric, radiometric or colorimetric detectors can be mounted onto the integrating sphere with an accompanying read-out instrument. In applications where wide variations in sphere housing temperature may occur, detectors with temperature stabilization are recommended to avoid drift of the measured signal. A spectrometer can also be used with or as an alternative to the integral light detector. Baffles are placed between the lamp and detector to prohibit direct illumination of



150 mm (6") dia top-load sphere with photometric detector for luminous flux measurement:

- 1. UMBB-150: 150 mm (6") dia Sphere
- 2. UMSS-SMT: sphere mount
- 3. UMSS-BT150: bench-top stand
- 4. UMPF-1.5: 1.5" dia port frame
- 5. UMLA-1.5/B: lamp adapter base
- 6. UMLA-150: lamp holder. Without socket
- UMPF-0.5: 0.5" dia port frame
 UMPA-0.5/11: port adapter for
- 11 type detector head

- 9. UMPA-0.5/RADIN: RADIN Diffuser for UMPA-0.5/11 port adapter
- 10. UMPB-PM: port baffle
- 11. TD-11VL01-2: temp. stabilized photometric detector
- 12. P-9710: single channel optometer with RS232 Interface
- 13. BN-0104: luminous flux transfer standard with UMLA-150 lamp holder
- 14. LPS-250: lamp power supply <u>Option:</u>
- P-2000: two channel benchtop optometer with RS232 and IEEE488 interface.

500 mm (20") dia sphere with photometric and UV-A detector for luminous flux and UV-A radiant power measurements.

- 1. UMBB-500: 500 mm (20") dia sphere
- 2. UMSS-HF500: hinge frame bench-top stand
- 3. UMPF-1.5: 1.5" dia port frame
- 4. UMLA-1.5/B: lamp adapter
- base 5. UMLA-500: lamp holder. Without socket
- 6. 2 x UMPF-0.5/11: 0.5" port frame
- 7. 2 x UMPA-0.5/11: port adapter for MD-11 detector head
- 8. 2 x UMPA-0.5/RADIN: diffuser

for UMPA-0.5/11 port adapter 9. UMPB-PM: port baffle

2

- 10. TD-11VL01-2: temp. stabilized photometric detector
- 11. TD-11UV01-2: temp. stabilized UV-A detector
- 12. P-2000: two channel benchtop optometer with RS232 and IEEE488 interface.
- 13. BN-0104: spectral radiant power transfer standard lamp 14. LPS-250: lamp power supply Options:
- UMLA-500: Additional lamp holder
- UMPA-2.0/AXL: Auxiliary lamp adapter



UM Series "Set-up Examples"

Uniform Light Source, Luminance and Spectral Radiance Standard

Uniform light sources are required to perform uniformity calibration of imaging systems and to calibrate luminance and spectral radiance measurement instruments.

Because of the integrating sphere's ability to produce multiple diffuse internal reflections and a uniform illumination of the inner sphere wall, it is the right tool to form the basis for uniform light sources.

Integrating spheres are available in different diameters offering the possibility for larger diameter emitting windows.

A number of optional accessory components can be added to the integrating sphere to expand its basic capabilities.

- Multiple lamps can be mounted onto one sphere offering high flexibility in intensity level.

If each lamps is operated with its own power supply, intensity can be varied in steps by simply switching light sources on and off.

- External spot lamps can be combined with variable attenuators to control the light output intensity in small steps. Attenuators with different resolutions are available

- External spot lamps can also be combined with manually exchangeable filters used for both neutral density attenuation and for forming the emission spectrum purposes.

- Remote control filter wheels and attenuators are another option to build a completely automated system.

- A large dynamic intensity range can be realized by combining neutral density filters with a variable attenuator

- Extra band pass filtering allows uniformity calibration across different wavelength ranges.

Large Aperture Uniform Light Source

- High speed shutters as fast as 40 ms enable calibration of the lamp's on/off characteristic.

- Luminance and spectral radi-

ance calibration is available through GO's calibration lab (ISO EN 17025 accredited for spectral irradiance & sensitivity).



Adjustable Uniform Light Source





Uniform light source 500 mm (20 in) dia sphere and (4) light sources for high uniformity and high luminance intensity. Color temperature controlled light sources

- 1. UMBB-500: 500 mm (20") dia. integrating sphere
- UMSS-SMT: sphere mount 3. UMSS-BT500: Bench top
- sphere stand
- UMPF-5.0: 5" dia port frame 5. UMPR-5.0/xx: port reducer. xx
- = 100 mm dia 6 4 x UMPB-SM: port baffle
- 4 x UMPF-1.5: 1.5" dia port 7
- frame

- 8. 4 x LS-IS1.5: light source with 50 W QH lamp
- 9. 4 x LPS-250: lamp power supply
- 10. 4 x BTH-19/2: 1/2 19" bench to housing for LPS-250
- 11. UMDP-37: detector port
- 12. UMPB-SM: port baffle
- 13. CT-3701: temperature stabi-
- lized color detector 14. P-9801: 8-channel meter
- 15. KLW: calibration of color temperature und luminance Options:

Software for remote control of intensity and color temperature control of each lamp

Adjustable uniform light source 300 mm (12 in) dia sphere and light source with manually controlled attenuator. Calibrated luminance detector for luminance intensity control

- 1. UMBB-300: 300 mm (12") dia integrating sphere
 - 2. UMSS-SMT: sphere mount
- 3. UMSS-300: bench top sphere stand
- 4. UMPF-2.0/3.5: port frame
- 5. UMPR-2.0/3.5/xx: port reducer, xx = 50 mm dia
- 6. UMPF-LSOK30: port frame 7. UMPB-PM: port baffle

- 8. LS-OK30-MIR: micrometer drive variable attenuator
- 9 LS-OK30-H100: spot lamp with 100 W QH lamp
- 10. LPS-250: lamp power supply 11. BTH-19/2: 1/2 19" bench top
- housing for LPS-250
- 12. UMPF-0.5: port frame
- 13. UMPA-0,5/11: port adapter
- 14. TD-11VL01-2: temperature stabilized photometric detector head
- 15. P-9710-1: optometer
- 16. KLW: calibration of luminance intensity

UM Series "Set-up Examples"

Multi-port Configurations for Reflectance, Transmittance, Absorbance

Integrating spheres offer unique functions which makes them very useful in reflectance. transmittance and absorbance applications

The Lambertian emittance characteristic of uniform light sources incorporating integrating spheres provides a hemispherical diffuse and uniform illumination of test samples e.g. such as light panels used for backlighting of displays.

In this case the luminance contrast of the light panels can be measured using an additional luminance meter viewing through an additional port.

The large acceptance angle of integrating sphere based light detection systems is the main reason for using spheres in reflectance and transmittance applications

Any light that is reflected by a

diffuse material or passing through a diffuse material must be completely detected. This should be independent of the angle of reflected or transmitted light and also from any increase in spot size. Flat field detectors do not meet this requirement. The port frame will not interfere with the acceptance angle if it is designed with a knife-edge.

Besides acceptance angle, larger measurement apertures can be made, an other argument for the integrating sphere.

A precisely collimated light beam is a prerequisite for accurate reflectance and transmittance measurements. The collimated beam must not only precisely target the sample area but also not contain any diffraction rings which can produce an offset light level inside the sphere limiting

measurement resolution. Consider the UP series integrating spheres for reflection, transmission, absorbance and other applications where more precise port positioning is critical.



Transmittance of Diffuse Transmitting Samples

25



Combined reflectance transmittance system built with two 150 mm (6 in) dia. spheres and two detectors.

- Reflectance sphere:
- 1. UMBB-150: 150 mm (6") dia integrating sphere
- 2. UMSS-14/M6: post stand
- 3. UMPF-1.0: 1.0" dia port frame
- 4. UMPF-1.5: 1.5" dia port frame
- 5. UMPR-1.5/xx: 1.5" port redu-
- cer with xx = 25 mm dia and
- knife edge
- 6. UMPB-SM: port baffle
- 7. UMPF-0.5: 0.5" dia port frame 8. UMPA-0.5/11: port adapter
- 9. UMPA-0.5/RADIN: diffuser

- 10. PD-1101-2: detector head
- Transmittance sphere:
- 11. UMBB-150: 150 mm (6") dia integrating sphere
- 12. UMSS-14/M6: post stand
- 13.UMPF-1.5: 1.5" dia port frame 14. UMPR-1.5/xx: 1.5" port redu-
- cer with xx = 50 mm dia and knife edge
- 15. UMPB-SM: port baffle
- 16. UMPF-0.5: 0.5" port frame 17. UMPA-0.5/11: port adapter
- 18. UMPA-0.5/RADIN: diffuser
- 19. PD-1102-2: detector head
- Option:
- 20. P-2000: 2-channel optometer

Reflectance sphere 500 mm (20") dia, 25 deg illumination port and photometric detector. Without bench top stand.

3

- 1. UMBB-500: 500 mm (20") dia integrating sphere
- 2. UMSS-SMT: sphere mount 3. UMPF-2.5/3.5: 2.5/3.5" port
- frame
- 4. UMPR-2.5/3.5/xx: 2.5/3.5" port reducer with xx = 50 mm diaand knife edge
- 5. UMPB-SM: baffle

Δ

6. UMDP-37: 37mm dia detector port

7. VL-3701-2: photometric illuminance detector with 37type package

8

8. UMPF-2.0: 2.0" dia port frame

Option:

- 9. P-9710: Optometer
- 10. LS-CB50: collimated light source set-up for 25 mm beam diameter
- 11. LPS-250: lamp power supply 12. BTH-19/2: 1/2 19" bench top housing

Hollow spheres with diameters from 100 mm (4 in.) to 1700 mm (68 in.) and three different coatings are offered to build integrating spheres to your need. Port frames, sphere stand and application accessory is assembled by your selection.

Please contact your local representation or the Gigahertz-Optik head office for personal assistance to discuss the set-up of the sphere and get your offer for the complete sphere set-up.



O.E.M. inquiries with different diameters, custom color painting of the sphere housing etc. are welcome.

UMB: Basic Spheres









UM Basic Spheres

Barium sulfate is an economical coating offering a \sim 97 % (555nm) diffuse reflectance characteristic within the wavelength range from 300 nm to 2400 nm. For photometric applications the reflectance can be reduce to \sim 80 % to reduce substitution effects errors produced by high reflectance.

Besides the standard ODP97 coating with best reflectance characteristic a more durable coating is available with slightly less reflectance capability.

OP.DI.MA. is Gigahertz-Optik's own white plastic with an excellent diffuse reflectance characteristic and high durability. Its diffuse reflectance is close to

that of a perfect lambertian reflector ensuring a highly uniform light distribution across the inner walls of the integrating sphere. The spectral reflectance covers a wide spectral range from ultraviolet to IR, 250 nm to 2500 nm with > 95 % reflectance and a peak reflectance of > 98 % in the VIS/NIR spectral range.

Uncoated gold outperforms all other uncoated metallic reflectors at infrared wavelengths. Gold's high reflectivity in the red portion of the visible spectrum is apparent by it's yellow hue. With a reflectance of \sim 95 % in the wavelength range from 800 nm to 20 μ m, gold is the best choice for surface coating of integrating spheres used in IR applications. A proprietary interior surface preparation creates the diffuse reflectivity required for uniform light distribution inside the sphere .

Ordering Information	Ordering Information & typical Specifications							
Model	Coating / Reflectance %	Hollow Sphere Diameter mm / in.	Sphere Frame Diameter mm / in.	Hollow Sphere Surface cm ²	5% Area Port Size mm / in.			
UMBB-100	ODP97 / 97%	100 / 4	124/4.9	314	43 / 1.7			
UMBB-150	ODP97 / 97%	150/6	174/6.8	707	66 / 2.5			
UMBB-210	ODP97 / 97%	210/8.3	234/9.2	1385	93 / 3.6			
UMBB-300	ODP97 / 97%	300/12	324/12.7	2827	130/5.1			
UMBB-500	ODP97 / 97%	500/20	524/20.6	7584	220/8.6			
UMBB-1000	ODP97 / 97%	1000 / 39	tbc	31416	450/17.7			
UMBB-1700	ODP97 / 97%	1700 / 68	tbc	90792	760/29.9			
UMBK-190	ODM98 / 98%	190/~7.5	234/9.2	1134	84 / 3.3			
UMBK-250	ODM98 / 98%	250/10	324/12.7	1964	112/4.4			
UMBK-460	ODM98 / 98%	460 / 18	524/20.6	6648	205/8.0			
UMBG-100	ODG95 / 95%	100/4	124/4.9	314	43 / 1.7			
UMBG-150	ODG95 / 95%	150/6	174/6.8	707	66 / 2.5			
UMBG-300	ODG95 / 95%	300/12	324 / 12.7	2827	130/5.1			

UM Series Components

UMSS-SM14, UMSS-SMT & UMSS-SMV Sphere Mount Adapters



The UMSS-SM14 sphere mount is the right choice for small diameter spheres up to 150 mm (6 in). The 14 mm diameter socket with M6 tapped hole allows the spheres to be mounted to commercial optical bench posts. The UMSS-SM14 is specified for upright use or for hanging.



The UMSS-SMT sphere mount is required to assemble UMB series spheres to the UMSS-BT bench-top stands using the M6 horizontal side holes provided. The SMT is fixed to the orbit flange. Additional M6 holes offer additional sphere mounting options into end-user applications.



The UMSS-SMV sphere mount enables spheres with a vertical orbit orientation to be assembled onto post mounts. The SMV is fixed to the orbit flange. Vertical M6 holes are provided for this purpose.



UMSS-SM14



 Ordering Information & typical Specifications

 Model
 Description
 Mounting Holes a

 UMSS-SM14
 Post mount socket fix mounted to sphere
 M6

 UMSS-SMT
 Pair of flange mounted triangular sockets fix mounted to sphere
 M6

 UMSS-SMV
 Pair of flange mounted sockets for fix mounted to sphere
 M6

Only available with UMB basic sphere order

UMSS-SMT UMSS-BT Bench-top Sphere Stands



The UMSS-BT stands are available for UMB series spheres from 150 mm (6") to 1700 mm (68") diameters.

The stands are made out of aluminum for lightweight stability. A modular design allows freedom for custom design modification or complete application setups.

Ordering Information & typical Specifications							
Stand Model	Sphere Model	Optical Axis (mm)	Width (mm)	Length	Al profile		
UMSS-BT150	UMBB&UMBK-150	118	220	150	20 x 20		
UMSS-BT210	UMBB-210, UMBK-190	tbc	tbc	tbc	20 x 20		
UMSS-BT300	UMBB-300, UMBK-250, UMBG-300	203	370	300	20 x 20		
UMSS-BT500	UMBB-500, UMBK-460	315	593	500	30 x 30		
UMSS-BT1000	UMBB-1000	tbc	tbc	tbc	tbc		
UMSS-BT1700	UMBB-1700	tbc	tbc	tbc	tbc		
Only available wi	h LIMB basic sphere order						

Only available with UMB basic sphere order

Page 102 - Integrating Spheres

UM Series Components

UMSS-HF Bench-top Hinged Frame Sphere Stands





The UMSS-HF bench-top sphere stand is designed for those applications requiring open access to the center of the integrating sphere for insertion and removal of the test lamp or sample. This is a common requirement in luminous flux or radiant power measurements of large multidirectional lamps, absorbance and transmission measurements. The UMSS-HF stands are designed for rugged long-term use. Port frames and other components like detector ports, light guide adapters, auxiliary lamps & lamps holders can be assembled to both sphere halves.

Ordering Information & Specifications							
Model	Sphere Model	Height of optical axis	Width	Length	Al profile		
UMSS-HF300	UMBB-300	tbc	tbc	tbc	tbc		
UMSS-HF500	UMBB-500	345	610	500	20 x 40		
UMSS-HF1000	UMBB-1000	tbc	tbc	tbc	tbc		
UMSS-HF1700	UMBB-1700	tbc	tbc	tbc	tbc		

Only available with UMB basic sphere order

UMPF Port Frames



Port frames are permanently mounted to the sphere at specified positions. Frames are available with free apertures from 0.5 to 5 inch diameter. A machined V-groove around the frame's



outer edge allows assembly of port plugs, reducers, adapters and other accessory components onto the port. Gold coated port frames are offered for use with UMBG series spheres.



Ordering Information & Specifications					
Model BaSO4	Model Gold	Free aperture diameter (a)	Outer diameter (b)	Frame Height (c)	
UMPF-0.5	UMPF-0.5-G	0.5 in. / 12.7 mm	0.75 in. / 18.8 mm	0.1 in. / 2.7 mm	
UMPF-1.0	UMPF-1.0-G	1.0 in. / 25.4 mm	1.25 in. / 31.5 mm	0.1 in. / 2.7 mm	
UMPF-1.5	UMPF-1.5-G	1.5 in. / 38.1 mm	1.75 in. / 44.2 mm	0.1 in. / 2.7 mm	
UMPF-2.0	UMPF-2.0-G	2.0 in. / 50.8 mm	2.25 in. / 56.9 mm	0.1 in. / 2.7 mm	
UMPF-2.0/3.5	UMPF-2.0/3.5-G	2.0 in. / 50.8 mm	3.5 in. / 88.9 mm	0.1 in. / 2.7 mm	
UMPF-2.5	UMPF-2.5-G	2.5 in. / 63.5 mm	2.75 in. / 69.6 mm	0.1 in. / 2.7 mm	
UMPF-3.0	UMPF-3.0-G	3.0 in. / 76.2 mm	tbc	0.1 in. / 2.7 mm	
UMPF-4.0	UMPF-4.0-G	4.0 in. / 101.6 mm	4.25 in. / 107.7 mm	0.1 in. / 2.7 mm	
UMPF-5.0	UMPF-5.0-G	5.0 in. / 127 mm	5.25 in. / 133.1 mm	0.1 in. / 2.7 mm	
Only available with UMB basic sphere order					

UMPP Port Plugs



Port plugs are used with UMPF port frames to close sphere port openings while maintaining it's uniformly coated interior surface. Plugs are fixed onto the side Vgroove of the port frame with three M2.5 screws. The surface facing the sphere interior is coated with barium sulfate or gold as applicable. Port plugs for UMBK spheres are available by special request.



Ordering Information & Specifications					
Model BaSO4	Model Gold 1)	Port Frame	Outer diameter a (mm / inch)	Height b (mm / inch)	
UMPP-0.5	UMPP-0.5-G	UMPF-0.5	27 / 1.06	5.3 / 0.2	
UMPP-1.0	UMPP-1.0-G	UMPF-1.0	39.9 / 1.57	5.3 / 0.2	
UMPP-1.5	UMPP-1.5-G	UMPF-1.5	52.6 / 2.07	5.3 / 0.2	
UMPP-2.0	UMPP-2.0-G	UMPF-2.0	65.3 / 2.57	5.3 / 0.2	
UMPP-2.0/3.5	UMPP-2.0/3.5-G	UMPF-2.0	97 / 3.81	5.3 / 0.2	
UMPP-2.5	UMPP-2.5-G	UMPF-2.5	78 / 3.07	5.3 / 0.2	
UMPP-3.0	UMPP-3.0-G	UMPF-3.0	tbc	5.85 / 0.23	
UMPP-4.0	UMPP-4.0-G	UMPF-4.0	116.1 / 4.57	5.85 / 0.23	
UMPP-5.0	UMPP-5.0-G	UMPF-5.0	141.5 / 5.57	5.85 / 0.23	
1) Only available with sohere system order					

UMPR Port Reducer



Port reducers for the UMPF series spheres enable the free aperture diameter of port frames to be variably decreased according to the application. The aperture diameter can be end-user specified up to the maximum aperture on order. The edge of the aperture is machined in a knife-edge design. Port reducers are fixed onto the V-groove of the port frame with three M2.5 screws. The surface facing the sphere interior is coated with barium sulfate or gold as applicable. Port plugs for UMBK



spheres are available by special request. Other optional forms of non-circular reducer such as slits can be custom ordered.

Ordering Information & Specifications								
Model BaSO4	Model Gold	Port Frame	Outer diameter a (mm / inch)	Height b (mm / inch)	xx = 0 to max free aperture *)			
UMPR-0.5/xx	UMPR-0.5/xx-G	UMPF-0.5	27 / 1.06	5.3 / 0.2	12.5 / 0.5			
UMPR-1.0/xx	UMPR-1.0/xx-G	UMPF-1.0	39.9 / 1.57	5.3 / 0.2	25.4 / 1.0			
UMPR-1.5/xx	UMPR-1.5/xx-G	UMPF-1.5	52.6 / 2.07	5.3 / 0.2	38.1 / 1.5			
UMPR-2.0/xx	UMPR-2.0/xx-G	UMPF-2.0	65.3 / 2.57	5.3 / 0.2	50.8 / 2.0			
UMPR-2.0/3.5/xx	UMPR-2.0/3.5/xx-G	UMPF-2.0	97 / 3.81	5.3 / 0.2	50.8 / 2.0			
UMPR-2.5/xx	UMPR-2.5/xx-G	UMPF-2.5	78 / 3.07	5.3 / 0.2	63.5 / 2.5			
UMPR-3.0/xx	UMPR-3.0/xx-G	UMPF-3.0	tbc	5.85 / 0.23	76.2 / 3.0			
UMPR-4.0/xx	UMPR-4.0/xx-G	UMPF-4.0	116.1 / 4.57	5.85 / 0.23	101.6 / 4.0			
UMPR-5.0/xx	UMPR-5.0/xx-G	UMPF-5.0	141.5 / 5.57	5.85 / 0.23	127 / 5.0			
*) free aperture dia	*) free aparture diameter win mm / inch with ar without knife wordre peeds to be apacified an order							

) free aperture diameter xx in mm / inch with or without knife wedge needs to be specified on of

UMDP Detector Ports



Detector ports are required to mount 37 and 45 mm diameter photometric, radiometric and color detectors to UM series integrating spheres.

The plastic ring that holds the detector electrically and thermally isolates the detector from the sphere.

It's free aperture can be in-



These port frames are needed to mount LS-OK30 and LS-OK34 external light sources to the integrating sphere. Typically these light sources are used together with the LS-OK30-VA variable attenuator.

creased according to the detector diffuser size. The sphere surface side of the detector adapters are coated with barium sulfate.

Simple adapters that fit into the detector port allow detector heads with different outside diameters to be mounted.



Ordering Information & Specifications						
Model	Detector dia.	Free aperture diameter *) a	Outside diameter b	Height c		
UMDP-37	37 mm	10 to 38 mm	46 mm	13 mm		
UMDP-45	42 & 45 mm	10 to 47 mm	56.9 mm	13 mm		

*) machined on order. Only available in order with UMB basic sphere

0 С

UMPF-LSOK Port Frames for Light Source LS-OK30 & LS-OK45

Ordering Information & Specifications					
Model	Aperture Dia. a	Dimension (mm)			
UMPF-LSOK30	30 mm	60 x 60			
UMPF-LSOK45	45 mm	tbc			
Only available with LIMP basic sphere order					

Only available with UMB basic sphere order

UMPB Port Baffles



Baffles are employed to shield ports from direct irradiation from an other port, a lamp, probe or from first reflection off the sphere surface.

UMDB-SM baffles are direct mounted to the sphere housing. They are typically used between two ports or between inside assembled lamps and ports.

UMDB-PM baffles are used in front of detector ports to prohibit direct irradiation of the detector from centrally mounted lamps. This baffles are typically post mounted.

UMDB-IL baffles are uses between to ports in-line assembled to the sphere. This arrangement is typically used for uniform light source set-ups with an external lamp at one port and the light

output at the other port.



UMPB-SM





UMPB-PM

UMPB-IL

Ordering Information				
Model BaSO4	Model Gold	Description		
UMPB-SM	UMPB-SM-G	Sphere mount baffle		
UMPB-PM	UMPB-PM-G	Post mount baffle		
UMPB-IL	UMPB-IL-G	In-Line spider mount baffle		
Only available with UMB basic sphere order				



Port adapters are universal tools used to assemble detectors. light guides, accessories and samples to UMPF sphere series port frames. The port adapter

UMPA -0.5/11

Port adapter specially designed to mount standard Gigahertz-Optik 11-type housing detectors or UFC fiber connectors with 11type adapter flanges to the UMPF-05 port frame.



can be fixed in place by three M2.5 screws to the port frame's standard V-groove. The port adapters can be replaced by port plugs to close the port as

UMPA -0.5/F & -1.0/F

Port adapter for use on UMPF-0.5 and UMPF-1.0 port frames respectively. The standard 1 mm dia. hole can be easily modified to larger diameters up to 0.5 / 1 in. (12.5 / 25.4 mm) to mount light guides and other accesso-



required by the application. The following list of port adapters represents only a small sample of available parts in Gigahertz-Optik's continuously expanding

program. Please contact your local representative or the factory for detailed drawings and also to keep abreast of new product developments.

UMPA Port Adapter

UMPA-4.0/QD1.9 & UMPA-4.0/QD3.0

Port adapter designed for uniform light source set-ups requiring a diffuser window. The 48 mm / 1.9 in dia. (model 1.9/QD)



or 76.2 mm / 3.0 in dia. (model 3.0/QD) diffuse quartz windows offer a useful wavelength range from 250-2500 nm.



UMPA-4.0/QD1.9

UMPA-4.0/QD3.0

UMPA-2.0/LED5 & UMPA-5.0/LED5

Port adapter which allows the precise positioning of LEDs mounted in the LED-5xx mount series of Instrument Systems /

Munich to the entrance port of Gigahertz-Optik's integrating spheres.



F





UMPA-2.0/AXL

Port adapter to set-up an auxiliary lamp by using a LS-IS1.5-H10 external light source. The adapter includes a hemisphere



baffle to protect ports and other components for direct illumination. The LS source can be direct





UMPA-2.0/AXL

LS-IS1.5 optional



UMPA-5.0/SH, SH-P & SH-CM

UMPA-1.5/LGA13

Port adapter to mount light guides with 13.2 mm dia. type adapter to the UM series integrating spheres.



5.0/SH-P port plug can also be used for custom modifications. The UMPA-5.0/SH-CM adapter is used when samples should be placed in the sphere center.

8



UMPA-5.0/SH

_В

Page 106 - Integrating Spheres

UM Series Components

Ordering Information & Specifications						
Model	Port Frame	Mechanical Interface	Outside Diameter	Height		
UMPA-0.5/11	UMPF-0.5	11-type detectors & UFC fiber couplers	25.4 mm / 1.0 in.	8.5 mm / 0.34 in.		
UMPA-0.5/11-RD	-	RADIN diffuser to be glued into UMPA-0.5/11	12 mm / 0.47 in.	1 mm / 0.04 in.		
UMPA-0.5/F	UMPF-0.5	1 mm dia., machineable up to 12.5 mm	25.4 mm / 1.0 in.	15 mm / 0.59 in.		
UMPA-1.0/F	UMPF-1.0	1 mm dia., machineable up to 25.4 mm	38.1 mm / 1.5 in.	15 mm / 0.59 in.		
UMPA-1.5/LGA13	UMPF-1.5	Socket for 13.2 mm dia light guide adapter	52.6 mm / 2.1 in	10.4 mm / 0.41 in		
UMPA-2.0/AXL	UMPF-2.0	Auxiliary lamp adapter with baffle for LS-IS1.5-10	65.3 mm / 2.6 in	6.4 mm / 0.25 in		
UMPA-2.0/LED5	UMPF-2.0	Socket for Instrument Systems LED-5xx LED adapter	65.3 mm / 2.57 in	25.85 mm / 1.01 in		
UMPA-4.0/QD1.9	UMPF-4.0	48 mm / 1.9 in dia quartz diffuser	116.1 mm/ 4.57 in	3.7 mm / 0.15 in.		
UMPA-4.0/QQ3.0	UMPF-4.0	76.2 mm / 3.0 in. diameter quart diffuser	116.1 mm / 4.57 in	3.7 mm / 0.15 in		
UMPA-5.0/LED5	UMPF-5.0	Socket for Instrument Systems LED-5xx LED adapter	141.5 mm / 5.57 in	25,85 mm / 1.01 in		
UMPA-5.0/SH	UMPF-5.0	101.6 mm / 4.0 in free aperture sample holder base	141.5 mm / 5.57 in	4.25 mm / 0.17 in.		
UMPA-5.0/SH-P	UMPF-5.0	Universal sample holder plate	141.5 mm / 5.57 in	6 mm / 0.24 in.		
UMPA-5.0/SH-CM	UMPF-5.0	Universal center mount sample holder	141.5 mm / 5.57 in	6 mm / 0.24 in.		

UMLA Lamp Adapter for Center Mount



UMLA lamp adapters are designed to mount lamps in the center of the integrating sphere for luminous flux or radiant power measurements.

The **UMLA-1.5/B** base offers four banana sockets for current and voltage connections to the lamp sample. It's rugged connector socket permits mounting and connection of the UMLA-xx lamp adapters. A UMPF-1.5 port frame is needed to mount the UMLA.1.5/B adapter onto the sphere.

The UMLA-xx lamp holders are

available in different lengths for compatibility with different size spheres. It's connector positioned at the top end of the UMLA-xx allows simple connection and mounting to the UMLA-1.5/B base socket. Lamp socket is supplied and assembled by the end-user.

The **UMLA-xx** includes a lamp socket specified by the end-user as a special order. Consult factory for assistance.



🕃 Gigahertz-Optik

Ordering Information & Specifications Model Accessory Description UMLA-1.5/B UMPF-1.5 Base adapter for UMLA-xx lamp holder **UMLA-300** Lamp holder for 300 mm / 12" diameter spheres UMLA-1.5/B **UMLA-500** UMLA-1.5/B Lamp holder for 500 mm / 20" diameter spheres UMLA-xx UMLA-1.5/B Lamp holder with lamp socket. Contact factory with detailed specs. Lamp adapter and lamp holder for UMBB-1000 and UMBB-1600 available by request

30 mm