

# Raman Spectroscopy

## Reliable, Non-Destructive Measurement

Ocean Optics offers a complete range of fully integrated Raman systems for handheld, laboratory and educational applications. Most of our Raman systems include a high-sensitivity spectrometer, a 532 nm or 785 nm laser, operating software and sampling accessories for probe- or cuvette-based analysis. Systems are compatible and portable.


Raman analysis is nondestructive, requires very limited sample preparation, and allows for sample volumes in the microliter range. In fact, Raman techniques can be used to acquire data through vial walls, pill pack windows and bags -- packaging forms that are frequently used in the biomedical and pharmaceutical industries.

Our Raman spectrometers are especially useful for identifying spectral features, or fingerprints, that are often distinct and clearly separated.

We offer off-the-shelf and modular approaches, sampling accessories for a variety of environments, cuvettes for solution measurements and a full range of probes for measuring solids and liquids.



Model	PeakSeeker 785	PeakSeeker Pro 785E	PeakSeeker Pro 785	PeakSeeker Pro 532	PinPointer	QE65000 for Raman
Laser Wavelength	785 nm	785 nm	785 nm	532 nm	785 nm	532 or 785 nm
System	Diode Laser Spectrometer Probe Laptop PC Software	Diode Laser Spectrometer Probe Laptop PC Software	Diode Laser Spectrometer Probe Laptop PC Software	Diode Laser Spectrometer Probe Laptop PC Software	Diode Laser Spectrometer Probe Software Embedded Computer with LCD Display	Choose from Laser Probes Sample Holders
Applications	Teaching labs, general purpose, moderate resolution Raman applications	Lab, industrial (pharmaceuticals), forensics	Lab, industrial (pharmaceuticals), forensics	Lab, industrial, security	Handheld, field-portable for on-site QC, forensics and security/inspection	General-purpose lab setups (aqueous solutions, powders, gels, tablets and surface media)
Item Code:	RAM-PKR-785	RAM-PRO-785E	RAM-PRO-785	RAM-PRO-532	RAM-PINPTR-785	QE65000, QE65000-RAMAN



### Technical Tip

Raman spectroscopy offers a number of benefits for testing and characterization. Like regular IR spectroscopy, it is rapid and non-destructive. Raman can capture data from a sample contained in plastic or other materials that are optically transparent to the wavelengths of interest.

Unlike IR spectroscopy, which falls spectrally within the water window, Raman spectroscopy can be used to capture data on aqueous samples or samples with high moisture content. And, with the emergence of economical diode lasers in the NIR spectral region, Raman is a more attractive

option for range of applications than ever before.

So, how to select between a Raman system with a 785 nm laser and one with a 532 nm laser? The 785 nm version is designed to minimize the fluorescence signal, making it useful for chemical identification and fingerprinting. What's more, the 785 nm version generates well defined peaks for use in chemometric analysis and can provide semi-quantitative data. Other wavelength options are available.

For applications where C-OH structural information is important, the 532 nm version is your best option. Such measurements are typical of biological and pharmaceutical sample analysis, where researchers study characteristics of active ingredients, binders, fillers and excipients.

# Benchtop Raman Systems


## Meet the PeakSeeker Pro™ Family



The PeakSeeker Pro family of benchtop Raman systems is accurate, cost-effective and easy to use.

These remarkable instruments use TE-cooled, high-efficiency CCD detector arrays and can be configured with either a 532 nm or 785 nm wavelength excitation laser.

PeakSeeker Pro is extremely versatile for measuring multiple sample types and comes with sampling accessories, software and USB connectivity that enables true plug and play operation.



### Versatility

You can sample solids, liquids and powders – even when samples are inside transparent packaging and containers such as bags, bottles and vials.

### Sensitivity

The PeakSeeker Pro fully meets the stringent requirements of USP Monograph 1120 for resolution, sensitivity and stability. The detector array is cooled to -20 °C. Deep blocking laser rejection filters obstruct Rayleigh scatter and isolate Raman scatter for valuable molecular analysis.

### Portability

The spectrometer is portable and lightweight. Rugged construction and ergonomic design allow it to be carried just about anywhere and set up for immediate use.

### Simplicity

Each PeakSeeker Pro system includes a notebook computer pre-loaded with RSIQ™ software for true point-and-click operation. Sample spectra are acquired within seconds.

Specifications	RAM-PKR-785	RAM-PRO-785	RAM-PRO-785E	RAM-PRO-532
Laser wavelength:	785 nm	785 nm	785 nm	532 nm
Laser power:	5-300 mW	5-300 mW	5-300 mW	100 mW
Spectral range, Raman shift:	200-3000 cm-1	200-2000 cm-1	200-3900 cm-1	200-4500 cm-1
Resolution:	~12 cm-1	6 cm-1	10 cm-1	10 cm-1
Signal to Noise Ratio:	250:1 (at full signal)	1000:1 (at full signal)	1000:1 (at full signal)	1000:1 (at full signal)
Wavelength stability:	<1 cm-1 in 12 hours	<1 cm-1 in 12 hours	<1 cm-1 in 12 hours	<1 cm-1 in 12 hours
Photometric stability:	<4% in 12 hours	<1% in 12 hours	<1% in 12 hours	<1% in 12 hours
Optics:	0.22 NA, sample spot size 0.1-.03 mm	0.22 NA, sample spot size 0.1-.03 mm	0.22 NA, sample spot size 0.1-.03 mm	0.22 NA, sample spot size 0.1-.03 mm
Power input:	15 VDC from power adapter input = 11/240 VAC @ 50/60 Hz	15 VDC from power adapter input = 11/240 VAC @ 50/60 Hz	15 VDC from power adapter input = 11/240 VAC @ 50/60 Hz	15 VDC from power adapter input = 11/240 VAC @ 50/60 Hz
Operating temperature:	-25 °C-+45 °C	-25 °C-+45 °C	-25 °C-+45 °C	-25 °C-+45 °C
Dimensions (cm):	36 x 29 x 11	36 x 29 x 11	36 x 29 x 11	36 x 29 x 11
Weight (kg):	3.6	3.6	3.6	3.6